Operating Instructions





Mabtech ASTOR™ 2

ELISpot reader

© 2023 Mabtech

Table of Contents

Introduction
Intended use
Symbols
Labels
Contact information
Regulatory information4
Quality Management System (QMS)5
Safety6
Product description8
Overview
In this package 8
Cleaning/Decontamination9
Disposal9
Warranty9
Technical specifications10
Setup and installation11
How to use the Mabtech ASTOR [™] 2 reader14
Service, maintenance, and transport16
To replace the fuses17

Introduction

Intended use

Mabtech ASTOR[™] 2: Detection of secreting cells in ELISpot assay. For research use only.

Symbols

These symbols are shown on the reader:

Symbols	Description
CE	Complies with the directives listed under "Regulatory information"
S/N	Serial number
	Electronic waste, disposal according to WEEE
Ń	Warning! The operating instructions must be consulted in all cases where this symbol is marked.

Labels Machine label:

> Mabtech ASTOR™ 2 ELISpot reader P/N: 2001 100-240 VAC, 50-60 Hz, 0.5 A max, 100 W Fuse: F4A H, 250 VAC, 5x20 mm Inputs: USB 3.0 2023-XX-XX Made in Sweden S/N: #XX Revision: XXX Box 1233, SE-131 28 Nacka Strand, SWEDEN www.mabtech.com



Contact information Mabtech AB Augustendalstorget 9 SE-131 52 Nacka Strand Sweden

Customer service: Phone: +46 8 716 27 00 E-mail: mabtech@mabtech.com

Regulatory information

The Mabtech ASTOR™ 2 ELISpot reader complies with the following directives:

Directive

Applied standards

CE mark for the European Economic Area (EEA)

2014/30/EU – Electromagnetic Compatibility Directive (EMC)	 EN IEC 61326-1, edition 4:2021-06-04 & edition 3:2013-01-11 EN 55011, edition 5:2016-04-01+A1:2017-04-21 +A11:2020-03-06+A2:2021-04-09 EN IEC 61000-3-2, edition 5:2019-03-01+A1:2021-04-09 & edition 4:2014-08-22 EN 61000-3-3, edition 3:2013-08-23+A1:2019-08-02+A2:2021-11-05
2014/35/EU – Low Voltage Directive (LVD)	 EN 61010-1, edition 3:2010-10-15+ A1:2019-02-22 EN 62471, edition 1:2008-09-12
2011/65/EU – Restriction of Hazardous Substances Directive (RoHS)	• CENELEC EN IEC 63000, edition 1:2018-12-07
2006/42/EC – Machinery Directive (MD)	
North America	
FCC Part 2	• §2.906 SDoC:2019-10-01
FCC Part 15*	 §15.107:2019-10-01 §15.109:2019-10-01
ICES-003	• Issue 7: 2020-10-15
International Standards	
EMC	• IEC 61326-1, edition 3:2020-10-26 & edition 2:2012-07-10

EMC	 IEC 61326-1, edition 3:2020-10-26 & edition 2:2012-07-10 	
	 CISPR 11, edition 6:2015-06-09+A1:2016-06- 23+A2:2019-01-18 	
	 IEC 61000-3-2, edition 5:2018-01-26+A1:2020-07-14 & edition 4:2014-05-26 	
	 IEC 61000-3-3, edition 3:2013-05-14+A1:2017-05- 18+A2:2021-03-25 	
Safety	• IEC 61010-1, edition 3:2010-06-10+A1:2016-12-16	
	 IEC 62471, edition 1:2006-07-26 	

*This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following

measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Quality Management System (QMS)

Standard	Description
EN ISO 13485:2016, Edition 4:2016-03-02 +A11:2021-09-08	 Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 9001:2015, Edition 5:2015-09-15	Quality management systems – Requirements

Safety

Before you install or operate the system, you must read the safety information. There is a risk of user injury or damage to the system if the instructions in this manual are not followed.



Caution! Only use the microtiter plate types specified under "Technical specifications." Other plates can cause malfunction and electrical failure.

(!) (!)

Caution! Do not insert microtiter plates that contain liquid. Risk of short circuit in the reader.

Caution! Do not place objects near or on the reader. Vibrations from the reader can make objects fall.

Product description



Overview

Number	Description	Number	Description
1	Cover	4	Programming switch*
2	Power button	5	USB port
3	Door to plate loader	6	IEC inlet with fuse holder

*Programming switch should be set downwards

In this package

The following parts are included in the package:

- Mabtech ASTOR™ 2 ELISpot reader
- USB cable
- Power cable
- Spare fuse
- Operating Instructions (this manual)
- MAIPSWU adapter

Cleaning/Decontamination

Clean the outside of Mabtech ASTOR[™] 2 ELISpot reader with a damp cloth when needed.

Before you send the system to any third party, fill out a certificate of decontamination to ensure that the instrument does not contain any potential contaminants.



Please note that Mabtech cannot accept any instrument or parts that may be contaminated with biohazardous substances such as viable biological agents, harmful quantities of chemicals, or radioactive materials.

Disposal



By the end of the reader life, recycle at a designated local reception point for Waste Electrical and Electronic Equipment (WEEE). You must not discard the reader as household waste.

Warranty

Refer to the Mabtech General Terms & Conditions of Sale found on http://www.mabtech.com/terms-and-conditions or the contractual agreement between the seller and the buyer.

Caution! Failure to obey the instructions in this manual will void the warranty.



Caution! The safety of any system or equipment that incorporates the Mabtech ASTOR ELISpot[™] 2 reader is the responsibility of the person connecting the reader to that system or equipment.

Technical specifications

Power supply

AC input voltage	100 to 240 VAC
AC input current	Max 0.5 A @100–240 VAC
AC frequency	50 to 60 Hz
Power consumption	40 W max
Fuses	Quick-Acting F4A H, 250 V(H = High Breaking Capacity (Ceramic Tube))

Operation and Maintenance

Interfaces	1 USB 3.0 port
Dimensions (H × W × D)	315 mm × 430 mm × 480 mm
Weight	19.1 kg
Environmental specifications	For indoor use only, +5°C to +40°C (+41°F to +104°F) Altitude up to 2000 m Relative humidity range: 5–95% Pollution Degree 2
Mains supply tolerance (%) or absolute mains supply values (voltage fluctuations of the nominal voltage)	+/-10%
Transient overvoltages	Overvoltage Category II Surge Level 3 (CM 12 ohm) & 4 (DM 2 ohm): 2 kV– 1.2/50 μs
IP class	IP20

Key features

Read/calculation speed	2 min / plate
Analyze speed	3 min from system start-up to excel export
Resolution (H × W)	1200 × 1200 pixels
Approved 96-well plate types	ELISpot MSIPS45 (ELISpot, clear) MSIPS4W (ELISpot, white) M8IPS4510 (ELISpot, strip) MAIPSWU10 (ELISpot, white, w/o underdrain MAIPS4510 (ELISpot, clear, w/o underdrain)

The Mabtech $\mathsf{ASTOR}^{\texttt{M}}$ 2 ELISpot reader must be connected to a computer equipped with the Mabtech $\mathsf{Apex}^{\texttt{M}}\text{-software}.$

Setup and installation



Warning! Follow the instructions in this manual when you lift or move the reader. Risk of injury to person or damage to the system.

Warning! Do not use or keep flammable liquids near the system. Risk of fire.

Caution! The reader is heavy. Place it on a solid table. Use two people to lift or move the reader.

Caution! Keep other equipment more than 10 cm away from the reader. Risk of malfunction caused by overheating.

Note: Keep the packaging in case the system needs to be sent back to Mabtech for service.



 Lift the Mabtech ASTOR[™] 2 ELISpot reader out of the box by the handles.



 Place the Mabtech ASTOR[™] 2 ELISpot reader upright on a solid table. Remove the inner packaging. Make sure that you have received all components specified in the section "In this package" on page 8. 3. Connect the USB and power cables.



Number	Description	Number	Description
1	Power cord ASTOR	5	Mouse cord
2	USB cord	6	Power cord monitor
3	HDMI/ DisplayPort cable	7	Power cord computer
4	Keyboard cord		

Warning! Do not connect faulty equipment to the system. Risk of electrical shock and damage to the system.

(!)

Caution! Connect the Mabtech ASTOR[™] 2 ELISpot reader to a grounded outlet. The input voltage must be between 100 to 240 VAC, 50-60 Hz.

i

Note: In the event of power failure, the reader might have to be re-started.



4. Push the button to turn on the Mabtech ASTOR[™] 2 ELISpot reader.



How to use the Mabtech ASTOR™ 2 reader

 Start the computer and double-click on the program symbol. The login screen is shown.



2. Select a profile and enter your password. Click Login.



3. Double-click on the type of assay you wish to run.



4. Click **Eject**. The reader ejects the plate loader.

+ New Plate	Open	D Read		Recount	← Eject	th sa	ave
Unnamed ×							
✓ Information							
ELISpot	A1	A2	A3	A4	A5	A6	A7

- **5.** Put the microtiter plate in the plate loader.
- 6. In Mabtech Apex[™] click **Read**. The reader retracts the plate loader. The microtiter plate is read and the spots are counted.



7. Click **Save** and select a folder. Click **Save**.

+ med	New Plate		Read Read Recount	Eject	Save 📝 Ex	port EG3 Prefe	prences	
	n JSpot ② Export To Folder					A6		
	← → * ↑ L Organize ▼ N	→ This ew folder	PC > Desktop > ASTOR VALIDATION >			O O Se	arch ASTOR VALIDATION	2
tory ajsa (on: LED ut	 Quick access Desktop Downloads Documents Pictures ASTOR VALIE HRP MT12A: TNF tinv for 	ATI B PR	Name 3490-3pw prod 9200504 LT 3490-3pw prod 6 annas HRP HRP MT12A3 PROD5 200430 LT2 KAJSA MONOCYT TITRERING TNF tiny for linearity	Date modified 2020-05-08 14:03 2020-05-08 14:00 2020-05-12 13:48 2020-05-08 15:14 2020-05-11 15:22	Type File folder File folder File folder File folder File folder	Size		
te or /well	File name: Save as type: A Hide Folders		Ng QC 3220-3PW prod 10			Sa	ve Cancel	

Service, maintenance, and transport



Caution! Only use shielded signal cables that are shorter than 3 m to comply with applicable EMC requirements. If the USB cable is longer than 2m, the connection between the reader and the computer can stop working.



Caution! Use only the parts specified below. Using other parts will void the warranty.

Return the equipment in the same box that it arrived in for repair/service/transport. Environmental conditions for transport and storage should be according to IEC 60068-2 and IEC 60721-3:

- Public transportation: -40 to +70°C
- Storage locations: -25 to +55°C.

Only the following parts can be replaced by the user:

Part	Description
Power cable	 240 V/10 A (international) or 100 V/15 A (North America) IEC 60320-1 C13 connector
USB cable	 Shielded USB cable Type A 3.0 male – type B 3.0 male ≤ 2m
Fuses x 2	 Cartridge fuses SP 5×20 250V 4A IEC 60127-2 Quick-Acting F4A, H (H = High Breaking Capacity (Ceramic Tube))



Note: All components in Mabtech ASTOR[™] are compliant with RoHS(Restriction of Hazardous Substances Directive, 2011/65/EU).

To replace the fuses



Warning! Disconnect the power before replacing the fuses. Risk of electrical shock.



- **1.** Pull out the power cable from the IEC inlet.
- 2. Remove the fuse holder from the IEC inlet with a flat screwdriver or similar tool.
- **3.** Replace the fuses.
- **4.** Plug the fuse holder and the power cable back in.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.







levices nent sys

The products are for research use only.

MABTECH shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages there from.

Mabtech AB (Head Office) Sweden Tel: +46 8 716 27 00 mabtech@mabtech.com

Mabtech, Inc. USA Tel: +1 513 871-4500 mabtech.usa@mabtech.com

www.mabtech.com