

GYROZEN遠心機

| セレクションガイド

		小型遠心機				多目的高速遠心機				低速遠心機	
		mini	1524	1536	1730R	1248	1248R	1580	1580R	406	416
冷却機能付					●		●		●		
最大遠心力											
アングル ローター	Max. 2,000xg									●	
	Max. 2,700xg										●
	Max. 12,300xg	●									
	Max. 14,980xg										
	Max. 15,520xg					●					
	Max. 16,500xg						●				
	Max. 21,200xg		●	●				●			
	Max. 25,900xg								●		
	Max. 27,200xg				●						
	Max. 32,300xg										
Max. 54,111xg											
スウィング ローター	Max. 2,800xg										●
	Max. 3,100xg										
	Max. 3,600xg										
	Max. 4,410xg										
Max. 5,300xg					●	●	●	●			
サンプル容器／ローターの種類											
PCRチューブ クライオチューブ	PCR strip	○	○	○	○						
	0.2/0.5ml	○	○	○	○	○	○	○	○		
	1.5/2.0ml	○	○	○	○	●	●	●	●		○
	Micro filter tube		○	○	○	○	○	○	○		
	1.8ml cryotube	○									
遠心管 大容量チューブ	3/4/5ml			○	○	●	●	●	●	○	○
	7ml					○	○	○	○		○
	10/15ml					●	●	●	●	○	●
	15ml conical					●	●	●	●	○	●
	50ml					●	●	●	●		○
	50ml conical					●	●	●	●		○
	85ml					●	●	●	●		○
	100ml					○	○	○	○		○
	100ml conical glass										
	250ml					○	○	○	○		
	250ml conical					○	○	○	○		
	500ml							○	○		
	500ml conical							○	○		
750ml							○	○			
1,000ml											
プレート	PCR plate					○	○	○	○		○
	Micro assay plate					○	○	○	○		○
	Tissue culture plate					○	○	○	○		○
	Deep well plate					○	○	○	○		○
	Filter set plate										
シリンジ	10/20/60ml					○	○	○	○		
血液バッグ	Double 200/400ml										

○ アングルローターのみ ○ スウィングローターのみ ● アングルローター、スウィングローターの両方とも可能

※価格に消費税は含まれておりません。